



FACIAL CARE



HIFU SKIN LIFTING



BODY SLIMMING



BREAST ENHANCEMENT



TATTOO REMOVAL



Jesús Galindo y Villa 3855 Casa 19  
Col. Cumbres del Nilo  
Guadalajara, Jalisco.

Tel. 01 (33) 1524-8420  
Cel. (33) 1360-8458, 33 1765-3715  
equiposdebellezagdl@gmail.com  
adyrahan@hotmail.com

[www.equiposdebellezagdl.webs.com](http://www.equiposdebellezagdl.webs.com)





**Mod. 1**  
 Name: PDT machine  
 Function: Seven LED colours,  
 Skin rejuvenation, wrinkle  
 removal, whitening  
 Voltage: 110V/220V 50-60Hz  
 Frequency: 1-30Hz  
 Working life: ≥5000 Hours  
 G.W.:20kgs  
 Package:86\*49\*36cm



**Mod. 2**  
 Name: PDT machine  
 Function: Seven LED colours,  
 Skin rejuvenation, wrinkle  
 removal, whitening  
 Voltage: 110V/220V 50-60Hz  
 Frequency: 1-30Hz  
 Working life: ≥5000 Hours  
 G.W.:14.5kgs  
 Package:80\*67\*46cm



**Mod. 3**  
 Name: PDT machine  
 Function: Four LED colours,  
 Skin rejuvenation, wrinkle  
 removal  
 Voltage: 110V/220V 50-60Hz  
 Frequency: 1-30Hz  
 G.W.:6kgs  
 Package:59\*31\*41cm  
 output Power: 12V



**Mod. 4**  
 Name: Face Neck LED Mask  
 Function: Blue and Red LED,  
 Skin rejuvenation, wrinkle  
 removal  
 Voltage: 110V/220V 50-60Hz  
 Power: 28W  
 G.W.:6.5kg  
 Package:40\*40\*30cm



**Mod. 5**  
 Name: Facial neck LED Mask  
 Function: 7 colours, EMS, Skin  
 rejuvenation, wrinkle removal  
 Voltage: 110V/220V 50-60Hz  
 Power: 15W  
 G.W.: 2kg  
 Package:30\*27\*22cm



**Mod. 6**  
 Name: Facial LED Mask  
 Function: 7 colours, Skin  
 rejuvenation, wrinkle removal  
 Voltage: 110V/220V 50-60Hz  
 Power: 7W  
 G.W.:1.5kg  
 Package:29\*25\*16cm



**Mod. 7**  
 Name: 3D HIFU machine  
 Function: 7 cartridges:3.0mm,4.5mm,  
 6.0mm, 8.0mm, 10mm,13mm,16mm;  
 11lines  
 Voltage: 110V/220V 50-60Hz  
 Power: 520W  
 Weight: 10kg  
 Package:50\*35\*45cm



### Mod. 8

Name: Diode lipo laser slimming  
 Function: 12pcs big and 2pcs small laser pads, slimming  
 Single-laser output:40mw-100mw  
 Power:650W  
 Voltage: 110V/220V 50-60Hz  
 G.W.:26kgs



### Mod. 9

Name: Diode lipo laser slimming  
 Function: 12pcs big and 2pcs small laser pads, slimming  
 Single-laser output:40mw-100mw  
 Power:650W  
 Voltage: 110V/220V 50-60Hz  
 G.W.:16kgs



### Mod. 10

Function: Cavitation, RF, Diode lipo laser slimming  
 Voltage: 110V/220V 50-60Hz  
 Power: 480W  
 G.W.:12kgs  
 Package:54\*47\*45cm  
 Single-laser output:40mw-100mw



### Mod. 11

Function: 10 big pads and 2 small pads, slimming  
 Voltage: 110V/220V 50-60Hz  
 Power: 480W  
 G.W.:12kgs  
 Package:54\*47\*45cm  
 Single-laser output:40mw-100mw



### Mod. 12

Function: 12pcs big and 2pcs small laser pads, slimming  
 Single-laser output:40mw-100mw  
 Power:650W  
 Voltage: 110V/220V 50-60Hz  
 G.W.:21kgs  
 Package:62\*46\*43CM



### Mod. 13

Function: 1. Multipolar RF Head  
 2. 40KHz Cavitation Head  
 3. 2pcs Vacuum probes  
 4. 10pcs big and 4pcs small laser pads  
 Single-laser output:40mw-100mw  
 Voltage: 110V/220V 50-60Hz  
 Power: 500W  
 G.W.:21kgs  
 Package:53\*50\*41CM



### Mod. 14

Function: 1, Tripolar RF head  
 2, 40KHz Cavitation Head  
 3, Vacuum Bipolar Head  
 4, Fourplar RF Head for face  
 5, Sixpolar RF Head for body  
 6. 6pcs big and 2pcs small laser pads  
 Single-laser output:40mw-100mw  
 Voltage: 110V/220V 50-60Hz  
 Power: 500W  
 G.W.:13kgs  
 Package:40\*39\*48cm



### Mod.15

Function: 1, Tripolar RF head  
 2, 40KHz Cavitation Head  
 3, Vacuum Bipolar Head  
 4, Fourplar RF Head for face  
 5, Sixpolar RF Head for body  
 Voltage: 110V/220V 50-60Hz  
 Power: 300W  
 G.W.:10kgs  
 Package:40\*39\*48cm



**Mod. 16**

Function: Slimming, skin Care  
 1, Twelve polar RF head  
 2, 40KHz Cavitation Head  
 3, Vacuum Bipolar Head  
 4, Fivepolar RF Head  
 5, BIO probe for face.  
 6pcs big and 2pcs small laser pads  
 Single-laser output: 40mw-100mw  
 Voltage: 110V/220V 50-60Hz  
 Power: 500W



**Mod. 17**

Function: Slimming, skin rejuvenation  
 1, Twelve polar RF head  
 2, 40KHz Cavitation Head  
 3, Vacuum Bipolar Head  
 4, Fivepolar RF Head  
 5, BIO probe for face  
 Voltage: 110V/220V 50-60Hz  
 Power: 300W  
 G.W.: 11kgs



**Mod. 18**

Function: Slimming, skin rejuvenation  
 1, Tripolar RF head  
 2, 40KHz Cavitation Head  
 3, Vacuum Bipolar Head  
 4, Sixpolar RF Head  
 5, Fourpolar RF Head  
 Voltage: 110V/220V 50-60Hz  
 Power: 300W  
 G.W.: 10kgs



**Mod. 19**

Function: Slimming, skin rejuvenation  
 1, Tripolar RF head  
 2, 40KHz Cavitation Head  
 3, Sixpolar RF Head  
 Voltage: 110V/220V 50-60Hz  
 Power: 120W  
 G.W.: 6.5kgs  
 Package: 44\*44\*28cm



**Mod. 20**

Function: Body slimming  
 40Khz cavitation probe,  
 Power: 120W  
 G.W.: 4kgs  
 Package: 35\*39\*23cm



**Mod. 21**

Function: 1. 2 ultrasonic probes for facial care  
 2. 40Khz Cavitation Head for body slimming  
 Voltage: 110V/220V 50-60Hz  
 Power: 75W  
 Weight: 5kgs  
 Package: 46\*37\*19cm



**Mod. 22**

Function: 1. 2 ultrasonic probes for facial care  
 2. 40Khz Cavitation Head for body slimming  
 Voltage: 110V/220V 50-60Hz  
 Power: 75W  
 Weight: 6kgs  
 Package: 46\*37\*29cm



**Mod. 23**

Function: 1. Cavitation probe for slimming  
 2. RF probe for fat dissolve  
 Voltage: 110V/220V 50-60Hz  
 Power: 110W  
 G.W.: 9.5kg  
 Package: 54\*49\*38cm



**Mod. 24**

Function: 1, Cavitation probe for slimming  
 2. RF probe for fat dissolve  
 3, 6pcs oriented cracking head  
 Voltage: 110V/220V 50-60Hz  
 Power: 150W  
 G.W.:12kg  
 Package:54\*49\*38cm



**Mod. 25**

Function: 1, Cavitation probe for slimming  
 2. RF probe for fat dissolve  
 3, 6pcs oriented cracking head  
 4. Ultrasonic probe for facial  
 Voltage: 110V/220V 50-60Hz  
 Power: 150W  
 G.W.:12kg  
 Package:42\*40\*50cm



**Mod. 26**

Function: Slimming, facial  
 1, Multipolar RF head  
 2, 40KHz Cavitation Head  
 3, Vacuum Bipolar  
 4, Bipolar RF for facial  
 5, BIO probe face lifting  
 Power: 250W  
 G.W.:24kg  
 Package:51\*46\*49cm



**Mod. 27**

Function: Slimming, facial  
 1, Tripolar RF head  
 2, 40KHz Cavitation Head  
 3, Vacuum Bipolar Head  
 4, Fourplar RF Head for face  
 5, Sixpolar RF Head for body  
 Voltage: 110V/220V 50-60Hz  
 Power: 300W  
 G.W.:21kgs



**Mod. 28**

Function: 1. Multipolar RF Head  
 2. 40KHz Cavitation Head  
 3. 2pcs Vacuum probes  
 4. 10pcs big and 4pcs small laser pads  
 Single-laser output:40mw-100mw  
 Voltage: 110V/220V 50-60Hz  
 Power: 500W



**Mod. 29**

Function: 1, Cavitation probe for slimming  
 2. RF probe for fat dissolve  
 3, 6pcs oriented cracking head  
 Voltage: 110V/220V 50-60Hz  
 Power: 150W  
 G.W.:31kgs



**Mod. 30**

Function: Slimming, facial  
 1, Multipolar RF head  
 2, 40KHz Cavitation Head  
 3, Vacuum Bipolar  
 4, Bipolar RF for facial  
 5, BIO probe face lifting  
 Power: 250W  
 G.W.:60kg



**Mod. 31**

Function: Slimming, facial  
 1, Tripolar RF head  
 2, 40KHz Cavitation Head  
 3, Vacuum Bipolar Head  
 4, Fourplar RF Head for face  
 5, Sixpolar RF Head for body  
 Voltage: 110V/220V 50-60Hz  
 Power: 300W  
 G.W.:60kgs



**Mod.32**  
 Name: Thermage Machine  
 Function: wrinkle removal, whitening  
 skin rejuvenation, face lifting  
 Voltage: 110V/220V 50-60Hz  
 RF frequency: 5MHZ  
 RF type: Thermage RF  
 Power: 350W  
 G.W.:7kgs



**Mod. 33**  
 Name: Thermage Freeze Machine  
 Function: wrinkle removal, whitening  
 skin rejuvenation, face lifting  
 Voltage: 110V/220V 50-60Hz  
 RF frequency: 5MHZ  
 RF type: Thermage RF  
 Power: 500W  
 G.W.:9.5kgs



**Mod. 34**  
 Name: Thermage Freeze Machine  
 Function: wrinkle removal, whitening  
 skin rejuvenation, face lifting  
 Voltage: 110V/220V 50-60Hz  
 RF frequency: 5MHZ  
 RF type: Thermage RF  
 Power: 500W  
 G.W.:20kgs



**Mod. 35**  
 Name: IPL machine  
 Function: hair removal, skin rejuvenation  
 pigmentation therapy,vascular  
 therapy, remove wrinkles  
 Voltage: 110V/220V 50-60Hz  
 Input Power: 1000W  
 Output Power:1800W  
 G.W.:30kg  
 Package:64\*52\*62cm



**Mod. 36**  
 Name: OPT IPL machine  
 Function: Fast hair removal, skin rejuvenation  
 pigmentation therapy,vascular  
 therapy, remove wrinkles  
 Voltage: 110V/220V 50-60Hz  
 Input Power: 1000W  
 Output Power:200W  
 G.W.:30kg  
 Package:64\*52\*62cm



**Mod. 37**  
 Name: RF +IPL machine  
 Function: hair removal, skin rejuvenation  
 pigmentation therapy,vascular  
 therapy, remove wrinkles  
 Voltage: 110V/220V 50-60Hz  
 Input Power: 1000W  
 Output Power:1800W  
 G.W.:35kg  
 Package:120\*66\*59cm



**Mod. 38**  
 Name: OPT IPL machine  
 Function: Fast hair removal, skin rejuvenation  
 pigmentation therapy,vascular  
 therapy, remove wrinkles  
 Voltage: 110V/220V 50-60Hz  
 Input Power: 1000W  
 Output Power:200W  
 G.W.:60kg



**Mod. 39**  
 Name: Portable IPL Machine  
 Function:hair removal, skin rejuvenation  
 pigmentation therapy,vascular  
 therapy, remove wrinkles  
 Voltage: 110V/220V 50-60Hz  
 Power: 200W  
 G.W.:2kgs  
 Package:36\*24\*20cm



**Mod. 40**

Name: HIFU Machine  
 Cartridges: 5 Heads(1.5mm, 3.0mm, 4.5mm, 8mm, 13mm)  
 Voltage: 110V/220V 50-60Hz  
 Power: 1500W  
 G.W.:20kg  
 Package: 48\*56\*58cm



**Mod. 41**

Name: HIFU Machine  
 Cartridges: 5 Heads(1.5mm, 3.0mm, 4.5mm, 8mm, 13mm)  
 Voltage: 110V/220V 50-60Hz  
 Power: 1500W  
 G.W.:40kg



**Mod. 42**

Name: HIFU Machine  
 Function: Vaginal Tightening  
 Voltage: 110V/220V 50-60Hz  
 Power: 1500W  
 G.W.:20kg  
 Probes: (3.0mm,4.5mm)



**Mod. 44**

Name: HIFU Machine  
 Function: Vaginal Tightening  
 Voltage: 110V/220V 50-60Hz  
 Power: 1500W  
 G.W.:40kg  
 Probes (3.0mm,4.5mm)



**Mod. 43**

Name: Radar Line tightening machine  
 Function: Facing Lifting, tighten Skin, Remove Wrinkle  
 Probe: 3.0mm and 4.5mm  
 Voltage: 110V/220V 50-60Hz  
 Power: 500W  
 G.W.:8kgs



**Mod. 46**

Name: 808 diode laser hair removal machine  
 Function: hair removal, skin rejuvenation  
 Voltage: 110V/220V 50-60Hz  
 Power: 2000W  
 G.W.:38kgs  
 Package: 135\*72\*65cm  
 Spot size: 10\*22mm



**Mod. 45**

Name: liposonix machine  
 Function: Body slimming  
 \*Shots Per probe: 525 shots  
 \* Depth: 1.3cm, 0.8cm  
 \* Energy: 5-75J/CM2  
 Voltage: 110V/220V 50-60Hz  
 Input Power: 200W  
 Output power: 500W  
 G.W.:6.5kg  
 Meas:45\*44\*28cm



### Mod. 47

Name: Hydra facial machine  
 Function: RF head  
 Oxygen injecting gun  
 Bubble skin cleaning pen  
 Micro current therapy head  
 Ley skin-repairing treatment head  
 Ultrasonic head  
 LEDMask  
 Voltage: 110V/220V 50-60Hz



### Mod. 48

Name: liposonix machine  
 Function: hair removal, skin rejuvenation, pigmentation therapy, vascular therapy, remove wrinkles  
 \* Shot quantity: 525 shots  
 \* Probes : 1.3CM, 0.8CM  
 \* Energy: 5-75J/CM2  
 Voltage: 110V/220V 50-60Hz  
 Input Power: 200W  
 Output power: 500W



### Mod. 49

Name: Hydra facial machine  
 Function: RF head  
 Oxygen injecting gun  
 Bubble skin cleaning pen  
 Micro current therapy head  
 Ley skin-repairing treatment head  
 Ultrasonic head  
 Voltage: 110V/220V 50-60Hz



### Mod. 50

Name: OPT+RF+ ND YAG Laser  
 Function: 1.OPT-SHR hair removal  
 2.cooling RF  
 3.ND YAG laser Tatoo removal  
 Voltage: 110V/220V 50-60Hz  
 power: 2500W  
 G.W.:55kgs  
 Package: 132\*65\*57cm  
 Filter: 480nm 530nm 640nm



### Mod. 51

Name: Tattoo removal laser  
 Function: 1.tattoo removal  
 2.tattoo removal on lip line, eyebrow, eyelid, body  
 3. pigment deposit removal  
 4 age spot, flat birthmark  
 Voltage: 110V/220V 50-60Hz  
 Power: 2000W  
 G.W.:16kgs  
 Package:56\*50\*40CM



### Mod. 53

Name: 4 in 1 infrared Pressotherapy Machine  
 Function: 10 pairs of air bags, 10 pairs of EMS, far infrared,  
 Intensity of pressure air wave: 75 kp.  
 Voltage: 110V /220V  
 Power: 500W  
 G.W.:52KGS  
 Packing: 120\*66\*59cm





### Mod. 52

Name: hydro oxygen injection machine  
 Function: Multi-polar RF handle  
 Cold hammer  
 micro needle  
 Oxygen injection  
 Voltage: 110V/220V 50-60Hz  
 Power: 500W  
 G.W.:20kgs  
 Package:53\*50\*41CM



**Mod. 54**

Name: pressotherapy lymphatic drainage  
 Intensity of pressure air wave: 75 kp.  
 Voltage: 110v /220v Power: 120W  
 G.W.: 20kgs  
 Air bag: 10 pairs



**Mod. 55**

Name: 3 in 1 Pressotherapy Machine  
 Function: 10 pairs of presso, 6 pairs of EMS, far infrared  
 Intensity of pressure air wave: 75 kp.  
 Voltage: 110V /220V Power: 500W



**Mod. 56**

Name: infrared Pressotherapy Machine  
 Function: 10 pairs of presso, far infrared  
 Intensity of pressure air wave: 75 kp.  
 Voltage: 110V /220V  
 Power: 750W



**Mod. 57**

Name: pressotherapy lymphatic drainage  
 Intensity of pressure air wave: 75 kp.  
 Voltage: 110v /220v  
 Power: 120W



**Mod. 58**

Name: Infrared air pressure slimming machine  
 Function: Infrared +Air pressure  
 Voltage: 110V/220V 50-60Hz  
 Power: 380W  
 weight:10kgs



**Mod. 59**

Name: infrared Pressotherapy Machine  
 Function: 10 pairs of presso, far infrared  
 Intensity of pressure air wave: 75 kp.  
 Voltage: 110V /220V  
 Power: 750W



### Mod. 60

Name: Breast Enlargement Equipment  
 Function: Breast Massage & Enlarger,  
 Cupping, Photon, Lymph detox,  
 Electro stimulation  
 Voltage: 110V/220V 50-60Hz  
 Power: 120W  
 weight: 8kgs



### Mod. 61

Name: Breast Enlargement Equipment  
 Function: Breast Massage & Enlarger,  
 Cupping, Photon, Lymph detox,  
 Electro stimulation  
 Voltage: 110V/220V 50-60Hz  
 Power: 120W  
 weight: 8kgs



### Mod. 62

Name: Professional Ultrasonic Machine  
 Function: 1. Eye Ultrasonic Probe  
 2. Facial Ultrasonic Probe  
 3. Body Ultrasonic Probe  
 Voltage: 110V/220V 50-60Hz  
 Power: 50W  
 G.W.:5kg  
 Package:46\*37\*19cm



### Mod. 63

Name: Needle Free mesotherapy  
 Function: melody muscle, promote  
 metabolism, decrease wrinkles,  
 whitening skin  
 Voltage: 110V/220V 50-60Hz  
 Power: 25W  
 G.W.:3.5kg  
 Package:39\*31\*23cm



### Mod. 64

Name: Needle Free mesotherapy  
 Function: 1. Whitening  
 2. tighten loose skin ;  
 3. Remove wrinkles  
 4. Treatment of acne ;  
 5. Skin Lifting, Skin Rejuvenation  
 Voltage: 110V/220V 50-60Hz  
 Power:100W



### Mod. 65

Name: Needle Free mesotherapy  
 Function: 1. Ultrasonic Probe  
 2. Nutrition Injection Probe  
 3. Freezing Lock Probe  
 4. Eye Injection Probe  
 5. BIO Wrinkle Dispelling Pen  
 Voltage: 110V/220V 50-60Hz  
 Power:65W



### Mod. 66

Name: Electro Muscle Stimulator  
 Function: Tighten body / Lose  
 weight / Chest Care  
 Voltage: 110V/220V 50-60Hz  
 Power: 35W  
 Weight: 4.5kgs  
 Package: 46\*36\*18cm



### Mod. 67

Name: Electro Muscle Stimulator  
 Function: Tighten body / Lose  
 weight / Chest Care  
 Voltage: 110V/220V 50-60Hz  
 Power: 35W  
 Weight: 4.5kgs  
 Package: 46\*36\*18cm



**Mod. 68**

Function: Tighten body  
Lose weight, Infrared Electro  
muscle stimulation pads  
negative electrode bag  
Voltage: 110V/220V 50-60Hz  
Power: 350W  
Weight: 6.5kgs  
Package: 46\*36\*18cm



**Mod. 69**

Function: Tighten body Lose  
weight, Infrared Electro  
muscle stimulation pads  
negative electrode bag  
Voltage: 110V/220V 50-60Hz  
Power: 500W  
Weight: 20kgs  
Package: 120\*63\*58cm



**Mod. 70**

Function: electro muscle stimulator  
G5 Vibration Massage  
Body Slimming, Breast Strengthen,  
Tighten skin  
Voltage: 110V/220V 50-60Hz  
Power:110W  
GW.:22kg  
Meas:102\*61\*46cm




**Mod. 71**

Name: Monopolar RET RF machine  
Function: Skin rejuvenation,  
wrinkle removal, whitening, face lifting  
Voltage:110V / 220V  
Power:300W  
RF frequency:0.3-0.5MHz  
RF energy:1-20J/cm3  
G.W.:6.5kg



**Mod. 72**

Name: Monopolar RF machine  
Function: Skin rejuvenation, wrinkle  
removal, whitening, face lifting  
Voltage:110V / 220V  
Power:300W  
RF frequency:0.3-0.5MHz  
RF energy:5-50J/cm3  
G.W.:6.5kg  
Meas:45\*37\*19cm



**Mod. 73**

Name: Monopolar RF machine  
Function: Skin rejuvenation, wrinkle  
removal, whitening, face lifting  
Voltage:110V / 220V  
Power:300W  
RF frequency:0.3-0.5MHz  
RF energy:5-50J/cm3  
G.W.:6.5kg  
Meas:44\*43\*28cm



**Mod. 74**

Name: Portable RF Machine  
Function:1 rf probe for face  
1 rf probe for body  
removal wrinkle, skin tightening  
Voltage: 110V/220V 50-60Hz  
Power: 100W  
G.W.:4.6kg  
Package:42\*34\*19cm



**Mod. 75**

Name: Portable RF Machine  
Function:3-polar for face and body  
removal wrinkle, skin tightening  
Voltage: 110V/220V 50-60Hz  
Power: 100W  
G.W.:4.6kg  
Package:42\*34\*19cm



**Mod. 76**

Name: Portable RF Machine  
 Function: 4-polar for face and body removal wrinkle, skin tightening  
 Voltage: 110V/220V 50-60Hz  
 Power: 100W  
 G.W.: 4.6kg  
 Package: 42\*34\*19cm



**Mod. 77**

Name: Portable RF Machine  
 Function: 6-polar for body removal wrinkle, skin tightening  
 Voltage: 110V/220V 50-60Hz  
 Power: 100W  
 G.W.: 4.6kg  
 Package: 42\*34\*19cm



**Mod. 78**

Function: RF skin tightening, Needle Mesotherapy, removal wrinkle, skin tightening  
 Voltage: 110V/220V 50-60Hz  
 Power: 110W  
 G.W.: 3kg  
 Package: 39\*29\*23cm



**Mod. 79**

Function: BIO electrodes, Magic gloves, face lifting  
 Voltage: 110V / 220V  
 Power: 35W  
 weight: 6.5kg  
 Meas: 45\*35\*18cm



**Mod. 80**

Name: G5 Vibration Massage  
 Function: Message / Dissolve fat  
 Voltage: 110V / 220V  
 Power: 100W  
 weight: 15.5kgs  
 Package: 102\*44\*42cm



**Mod. 81**

Name: G5 Vibration Massage  
 Function: Message / Dissolve fat  
 Voltage: 110V / 220V  
 Power: 35W  
 Weight: 6kg  
 Package: 50\*35\*21cm



**Mod. 82**

Name: 8 in 1 multifunctional  
 Function: Ultrasonic / Brush Vacuum / Spray / Remove spot / Galvanic / High-frequency  
 Voltage: 110V/220V 50-60Hz  
 Power: 35W  
 Weight: 9.5kg  
 Package: 65\*37\*19cm



**Mod. 83**

Name: 8 in 1 multifunctional  
 Function: Ultrasonic / Brush Vacuum / Spray / Remove spot / Galvanic / High-frequency  
 Diamond microdermabrasion  
 Hot therapy / Cold therapy  
 Voltage: 110V/220V 50-60Hz  
 Power: 100W  
 G.W.: 13kgs  
 Package: 40\*39\*48cm



**Mod. 84**  
 Function: RF vacuum, Lymph drainage, cupping  
 Voltage: 110V/220V 50-60Hz  
 Power: 100W  
 G.W.: 6kgs  
 Package: 42\*34\*19cm



**Mod. 85**  
 Function: BIO electrodes  
 Magic gloves, Hot cold hammer  
 Voltage: 110V/220V 50-60Hz  
 Power: 30W  
 weight: 6.8kg  
 Package: 50\*35\*22cm



**Mod. 86**  
 Function: 2 zone Infrared thermal blanket, Reduce weight, Detox  
 Voltage: 110V/220V 50-60Hz  
 Power: 650W  
 Weight: 7.5kgs  
 Package: 82\*48\*16cm



**Mod. 87**  
 Function: 3 zone Infrared thermal blanket, Reduce weight, Detox  
 Voltage: 110V/220V 50-60Hz  
 Power: 850W  
 Packing : 62\*46\*43cm  
 weight: 12kg



**Mod. 88**  
 Function: Patter breast care, nipple care, breast enlargement  
 Voltage: 110V/220V 50-60Hz  
 Power: 70W  
 weight: 6.6kg  
 Package: 37\*19\*64cm



**Mod. 89**  
 Name: Electro Muscle Stimulator  
 Function: Tighten body / Lose weight  
 E-bags  
 Voltage: 110V/220V 50-60Hz  
 Power: 50W



**Mod. 90**  
 Function: Microdermabrasion, facial peeling, skin rejuvenation  
 Voltage: 110V/220V 50-60Hz  
 Power: 50W  
 weight: 5kg  
 Package: 37\*46\*15cm



**Mod. 91**  
 Function: Microdermabrasion, facial peeling, skin rejuvenation  
 Voltage: 110V/220V 50-60Hz  
 Power: 50W  
 weight: 7kg  
 Package: 37\*49\*19cm



### Mod. 92

Function: 1. Microdermabrasion  
 2. Ultrasonic  
 3. Skin scrubber  
 4 & 5. Hot & cool hammer  
 Voltage: AC 110V/220V 50-60Hz  
 Power: 110W  
 G.W.: 10.5kgs  
 Package: 48\*39\*31cm



### Mod. 93

Function: 1. 2 in 1 crystal powder  
 2. diamond Microdermabrasion  
 Voltage: AC 110V/220V 50-60Hz  
 Power: 150W  
 G.W.: 18.5KG  
 Package: 57\*46\*38cm



### Mod. 94

Function: Cryolipolysis machine  
 with 1 Handle  
 Voltage: AC 110V/220V 50-60Hz  
 Power: 500W  
 Pressure output: 0-100KPa  
 Cryo handle size : 140X50mm  
 G.W.: 28KGS  
 Package: 66\*48\*41cm



### Mod. 95

Function: 1. Cryolipolysis machine  
 2. Vacuum Roller  
 Voltage: AC 110V/220V 50-60Hz  
 Power: 500W  
 Pressure output: 0-100KPa  
 Cryo handle size : 140X50mm  
 G.W.: 28KGS  
 Package: 66\*48\*41cm



### Mod. 96

Function: 1. Cryolipolysis probe  
 2. Vacuum Roller  
 3. Lipo Laser  
 4. Cavitation  
 5. RF  
 Voltage: AC 110V/220V 50-60Hz  
 Power: 500W  
 Pressure output: 0-100KPa  
 Cryo handle size : 140X50mm  
 G.W.: 32KGS  
 Package: 66\*58\*51cm



### Mod. 97

Function: Two Cryolipolysis probes  
 Voltage: AC 110V/220V 50-60Hz  
 Power: 500W  
 Pressure output: 0-100KPa  
 Cryo handle size : 140X50mm  
 G.W.: 20KGS  
 Package: 66\*48\*41cm



### Mod. 98

Function: Ozone hot Steamer  
 Voltage: AC 110V/220V 50-60Hz  
 Power: 750W  
 Weight: 8.5kg  
 Package: 90\*43\*30cm



### Mod. 99

Function: Ozone Steamer and  
 magnify lamp  
 Voltage: AC 110V/220V 50-60Hz  
 Power: 750W  
 G.W.: 11.5kg  
 Package: 90\*43\*41cm



### Mod. 100

Function: High-frequency, Galvanic  
Voltage: 110V/220V 50-60Hz  
Power: 35W  
weight: 5KGS  
package: 50\*37\*19cm



### Mod. 101

Function: Remove spot  
Voltage: 110V/220V 50-60Hz  
Power: 12W  
G.W.: 3.5kg  
Package: 50\*37\*19cm



### Mod. 102

Function: Remove spot,  
Electrotherapy, Vacuum,  
Spray, Galvanic  
Voltage: 110V/220V 50-60Hz  
Power: 65W  
weight: 7kg  
package: 58\*37\*19cm



### Mod. 103

Function: Ultrasonic, Vacuum,  
Spray, High-frequency,  
Remove spot, Breast care  
Voltage: 110V/220V 50-60Hz  
Power: 120W  
weight: 8kg  
Package: 65\*37\*19cm



### Mod. 104

Function: Ultrasonic, Vacuum,  
Spray, High-frequency  
Voltage: 110V/220V 50-60Hz  
Power: 80W  
G.W.: 6kgs  
Package: 50\*37\*19cm



### Mod. 105

Function: Ultrasonic,  
BIO skin lifting  
Voltage: 110V/220V 50-60Hz  
Power: 50W  
weight: 6kgs  
package: 49\*37\*19cm



### Mod. 106

Function: Steamer/ Ultrasonic  
Ozone/Skin scrubber/ Cold &  
Hot therapy/ Brush/ High-frequency  
Galvanic/ Vacuum/ Spray/ Remove  
spot/ Magnifying Lamp/ Woods Lamp  
Voltage: 110V/220V 50-60Hz  
Power: 950W  
G.W.: 31kg  
Package: 110\*47\*54cm



### Mod. 107

Function: Hot ozone steamer/Ozone  
Ultrasonic/ Brush/ Highfrequency/  
Galvanic/ Vacuum/ Spray/ Remove  
spot/ Magnify lamp/ Woods lamp  
Voltage: 110V/220V 50-60Hz  
Power: 750W  
G.W.: 38.5kg  
Package: 104\*77\*51cm



**Mod. 108**

Function: two ultrasonic  
 Voltage: 110V/220V 50-60Hz  
 Power: 25W  
 weight: 1.6kg  
 Package:35\*11\*24cm



**Mod. 109**

Function:  
 Ultrasonic and Remove spot  
 Voltage: 110V/220V 50-60Hz  
 Power: 25W  
 weight: 1.6kg  
 Package:35\*11\*24cm



**Mod. 110**

Function: Boxy Skin &  
 Hair Analyzer  
 Voltage: 110V/220V 50-60Hz  
 G.W.:1.86kg  
 Package:33\*14\*26cm



**Mod. 111**

Function: Boxy Skin &  
 Hair Analyzer  
 Voltage: 110V/220V 50-60Hz  
 Power: 55W  
 G.W.:5kg  
 Package:38\*30\*26cm



**Mod. 112**

Function:  
 Hot Ozone Steamer  
 Ultrasonic  
 Skin scrubber  
 Cold & Hot Therapy  
 Brush  
 High-Frequency  
 Galvanic  
 Vacuum  
 Spray  
 Remove spot  
 Magnifying Lamp  
 Woods Lamp  
 Vacuum roller  
 Micro-dermabrasion  
 Voltage: 110V/220V 50-60Hz  
 Power: 950W  
 G.W.:31kg  
 Package:110\*47\*54cm



**Mod. 113**

Function:  
 Hot steamer  
 Ultrasonic  
 High-frequency  
 Magnifying Lamp  
 Voltage: 110V/220V 50-60Hz  
 Power: 750W  
 G.W.:23kg  
 Package:101\*50\*50cm



**Mod. 114**  
 micro needle pen  
 Package:21\*15\*9cm



**Mod. 115**  
 High frequency device  
 Package:34\*21\*6cm



**Mod. 116**  
 Diamond Peeling  
 Package:27\*21\*11cm



**Mod. 117**  
 Function: Freckle removal,  
 skin whitening,  
 Fantasy wrinkle removal  
 Voltage: AC 110V/220V 50-60Hz  
 Power: 25W  
 G.W.:0.5kgs  
 Package:20\*17\*6cm



**Mod. 118**  
 Function: Thermage wrinkle removal,  
 skin rejuvenation  
 Needle Length: 0.25/0.5/0.75/1.0/  
 1.5/2.0/2.5/3.0(mm)  
 G.W.: 0.5kgs  
 Package: 27\*20\*9cm